Sterling Series  
**Sprague  
Rappaport-Type  
Stethoscope**

**INSTRUCTION MANUAL**

English • Spanish

**Latex-Free**

---

# Table of Contents

---

Introduction .....	4
Product Features .....	4
Ensuring a Proper Fit.....	5
Replacing Parts .....	6-7
Cleaning .....	8
Care and Safety Information .....	8
Instrucciones en Español.....	9

---

**Toll-Free Customer Care Help Line: 1-866-326-1313**  
**Monday – Friday 8:30 a.m. – 4:30 p.m. CST**


Distributed by  
Veridian Healthcare  
1175 Lakeside Drive  
Gurnee, IL 60031  
[www.veridianhealthcare.com](http://www.veridianhealthcare.com)  
Made in China  
#93-1049 03/15  
©2015 Veridian Healthcare

---

## Product Features

---

---

## Introduction

---

Thank you for purchasing this Veridian Sterling Series Sprague Rappaport-Type Stethoscope. Please read this instruction manual completely, to ensure safe use, accurate results and long life of the stethoscope.

This stethoscope is designed to deliver the most reliable and precise performance of any stethoscope on the market. Its ease of use and quality design and construction are uncompromising; this is a top-tier tool for cardiologists and other medical professionals.

With proper care, your new stethoscope will provide you and your patients with many years of safe, accurate and comfortable performance.

---

---

## Ensuring a Proper Fit

---


The earpieces on your Veridian Healthcare Stethoscope have been carefully calibrated to provide clear and audible acoustics. To best hear the sounds to be analyzed, please keep the following tips in mind:

Multiple earpieces are included. Select the pair most comfortable and effective for you. The earpieces should point forward when in use, *fig. 1*.

The earpieces should fit in your ears securely but comfortably. The binaurals may be adjusted to ensure maximum comfort and efficacy. To increase the tension, simply squeeze the binaurals just above the joint where they meet, crossing the earpieces, *fig 2*. Do not apply too much pressure; apply a minimum of force, slightly increasing the pressure until the earpieces fit comfortably.


Should the earpieces be too tight, pull the binaurals apart, grasping them just above the joint where they meet, *fig. 3*. Be gentle, so as not to damage the stethoscope.

Always use caution when adjusting the binaurals to avoid damage. When bringing the earpieces closer together, the maximum tension used should bring the black tubing on either side within contact. When pulling the binaurals apart, do not exceed 90 degrees between the binaurals.


*Fig. 1*  
Correct  
Position

*Fig. 1*  
Incorrect  
Position

*Fig. 2*

*Fig. 3*

# Replacing Parts

Note: Use only Veridian replacement parts in this unit, available from your retailer or directly from Veridian Healthcare at 866-326-1313.

## Diaphragm and Bell Use:



*Large Diaphragm* - Detects diastolic, first, second and third heart sounds, low frequency heart sounds and high-pitched murmurs in adults.

*Small Diaphragm* - For high-pitched infant or pediatric heart sounds.

*Large Bell* - Pinpoints low and medium pitched sounds and murmurs in adults.

*Medium Bell* - Detects low and medium pitched sounds or murmurs in narrow or small places.

*Small Bell* - For low and medium pitched sounds and murmurs from an infant or pediatric patient.

---

# Replacing Parts

---

## **Changing the Diaphragm:**

Veridian Healthcare Sprague Rappaport-Type Stethoscopes include a chestpiece key for easy diaphragm replacement, *fig 1*.


Remove either the small or large diaphragm from the chestpiece.

Place the chestpiece diaphragm side down on the palm of your hand. Using the key, twist and remove the inner metal disc for access to the plastic diaphragm, *fig 2*.


After replacing the diaphragm, use the metal key to tighten the metal disc into place and reattach to the chestpiece, *fig 3*.

## **Changing the Components:**


Anytime a component is changed, tap one side of the chestpiece to confirm that the stethoscope is functioning properly. If no sound is heard through the binaural tips, rotate the head 180° or until a clicking sound is heard. Repeat tap.

*Fig. 1*

*Fig. 2*

*Fig. 3*

---

## Cleaning

---

The entire unit may be cleaned with either soap and water or isopropyl alcohol. The earpieces may be removed from the binaurals to facilitate cleaning. After each cleaning, all parts should be thoroughly rinsed and dried completely. Never immerse any part of the stethoscope in liquid. Never use solvents, oils or sterilization equipment to clean the unit.

The chestpiece may be cleaned with the included cleaning cloth as needed. Soapy water or isopropyl alcohol may also be used.

The diaphragm may also be removed and cleaned with soapy water or isopropyl alcohol. Rinse and dry completely after cleaning.

---

## Care and Safety Information

---

Improper care or use of your stethoscope may result in damage to the unit. Following these instructions will ensure the stethoscope's accuracy and long life.

- Do not store the stethoscope in direct sunlight, dusty or humid environments, or extreme temperatures.
- Never use solvents or abrasive cleaning products to clean the stethoscope. Do not immerse the stethoscope or use oils to clean it. Follow the cleaning instructions provided below.

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