Table of Contents

Introduction	3
Product Features	4-5
TrueTone Technology.....	6
Ensuring a Proper Fit.....	7
Replacing Parts	8
Cleaning	9
Care and Safety Information	9
Instrucciones en Español.....	10-18

Introduction

Thank you for purchasing a Veridian Pinnacle Series Stethoscope. Please read this instruction manual completely, to ensure safe use, accurate results and long life of the stethoscope.


This stethoscope is designed to deliver the most reliable and precise performance of any stethoscope on the market. Its ease of use and quality design and construction are uncompromising; this is a top-tier tool for cardiologists and other medical professionals.

With proper care, your new stethoscope will provide you and your patients with many years of safe, accurate and comfortable performance.


Toll-Free Customer Care Help Line: 1-866-326-1313
Monday – Friday 8:30 a.m. – 4:30 p.m. CST

Product Features - Low Profile Stethoscopes

Product Features - Dual Head Stethoscopes

TrueTone™ Technology - Low Profile ONLY

Your stethoscope uses an advanced TrueTone diaphragm, setting it apart from previous generations of stethoscopes. TrueTone technology allows both high- and low-frequency sounds to be clearly heard without changing or repositioning the chestpiece.

To hear low-frequency sounds, simply place the chestpiece firmly against the patient's skin, using light pressure, *fig 1*.

To listen for high-frequency sounds, apply more pressure, pressing the chestpiece firmly against the patient's chest, *fig. 2*. The diaphragm will adjust to pick up these higher-frequency sounds.


Fig. 1


Fig. 2

Ensuring a Proper Fit

The earpieces on your Veridian Healthcare Stethoscope have been carefully calibrated to provide clear and audible acoustics. To best hear the sounds to be analyzed, please keep the following tips in mind:

Multiple earpieces are included. Select the pair most comfortable and effective for you. The earpieces should point forward when in use, *fig. 1*.

The earpieces should fit in your ears securely but comfortably. The binaurals may be adjusted to ensure maximum comfort and efficacy. To increase the tension, simply squeeze the binaurals just above the joint where they meet, crossing the earpieces, *fig 2*. Do not apply too much pressure; apply a minimum of force, slightly increasing the pressure until the earpieces fit comfortably.

Should the earpieces be too tight, pull the binaurals apart, grasping them just above the joint where they meet, *fig. 3*. Be gentle, so as not to damage the stethoscope.

Always use caution when adjusting the binaurals to avoid damage. When bringing the earpieces closer together, the maximum tension used should bring the black tubing on either side within contact. When pulling the binaurals apart, do not exceed 90 degrees between the binaurals.


Fig. 1
Correct
Position

Fig. 1
Incorrect
Position



Fig. 2

Replacing Parts

Note: Use only Veridian replacement parts in this unit, available from your retailer or directly from Veridian Healthcare at 866-326-1313.

Replacing the Diaphragm:

1. Pull the plastic retaining ring, *fig. 1* off of the top, metal portion of the chestpiece. Take note that the diaphragm side of the ring is round, while the chestpiece side is flat.
2. Gently remove the diaphragm from within the ring.
3. Insert the new diaphragm into the ring, through the rounded (diaphragm side) edge. Ensure that the diaphragm is placed securely in the groove of the ring.
4. Roll the ring, with the diaphragm securely in place, into the chestpiece, again ensuring a snug fit. The ring should be securely in place, yet allow easy rotation. The ring should not be easily separated from the chestpiece.
5. Check once again that the diaphragm remains secure and snug within the ring, and the ring is firmly in place in the chestpiece. Should there be any problems, make adjustments or start over.


Fig. 1

The earpieces may also be replaced, as needed. To replace the earpieces, unscrew them from the binaurals and screw the replacement earpieces on securely.

Cleaning

The entire unit may be cleaned with either soap and water or isopropyl alcohol. The earpieces may be removed from the binaurals to facilitate cleaning. After each cleaning, all parts should be thoroughly rinsed and dried completely. Never immerse any part of the stethoscope in liquid. Never use solvents, oils or sterilization equipment to clean the unit.

The chestpiece may be cleaned with the included cleaning cloth as needed. Soapy water or isopropyl alcohol may also be used.

The diaphragm may also be removed and cleaned with soapy water or isopropyl alcohol. Rinse and dry completely after cleaning.

Care and Safety Information

Improper care or use of your stethoscope may result in damage to the unit. Following these instructions will ensure the stethoscope's accuracy and long life.

- Do not store the stethoscope in direct sunlight, dusty or humid environments, or extreme temperatures.
- Never use solvents or abrasive cleaning products to clean the stethoscope. Do not immerse the stethoscope or use oils to clean it. Follow the cleaning instructions provided below.

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