

### **Product Information**

## ISO Class 5 - 8

Cleanroom Class 100 – 100,000 EU Grade A-D USP <797>

# SatPax® 1000

## Pre-wetted 55% Cellulose / 45% Polyester Nonwoven Cleanroom Wiper

SatPax® 1000 combines Durx® 670 nonwoven wipers composed of a hydroentangled nonwoven blend of 55% cellulose and 45% polyester with various saturation levels of isopropyl alcohol and DI water.

This pre-wetted format provides a cost effective and easy to use solution versus traditional bulk handling of solvents, maintenance of squirt bottles and inconsistent wetting and cleaning associated with wetting a dry wiper.



Other Class 5 and above Pre-wetted wipers

SatPax® 670-R SatPax® 550

#### Kev Attributes

- 55% cellulose / 45% polyester hydroentangled nonwoven blend
- No chemical binders in base material
- Pre-wetted with consistent IPA / DI Water concentrations and saturation levels
- Re-sealable solvent resistant packaging with wipers in a c-folded configuration for single withdrawal
- Also available in canister dispensers—wipers in perforated roll configuration for pop-up convenience

#### **Benefits**

- Low extractables and fiber and particle counts
- Smooth and durable with good wet strength
- Reduces alcohol usage and preparation / handling costs
- Reduces VOC emissions
- Increases cleaning efficiency
- Increases cleaning protocol consistency

#### Environmenta

- Reduces VOC (Volatile Organic Compound) emissions
- Complies with RoHS and pre-registered under REACH

#### Applications

- Designed for use in ISO Class 5 and higher cleanroom environments and USP <797> applications
- Designed for use in wet cleaning of critical surfaces where control of flammable solvents and flammable solvent concentrations is required
- Final cleaning of surfaces or products prior to manufacturing or packaging
- High saturation level is ideal for removing cleaning and disinfecting residues in regulated environments

#### **Alcohol Mixtures**

Alcohol / DI Water mixtures can be varied to fit the customer requirements. Typical mixtures are 70/30 and 9/91 IPA/DI Water.

#### **Saturation Levels**

The amount of solution contained in each wiper can be varied according to customer requirements. Higher saturation levels apply more solution to the surface during cleaning.

#### Sterile Validated Option

For aseptic processing areas, the same wiper material can be provided in a gamma irradiated validated sterile to a 10<sup>-6</sup> sterility assurance level.

Americas	Tel 1 413 528 2602	Fax 1 413 528 2614	usa@berkshire.com
Europe	Tel 44 (0) 870 757 2877	Fax 44 (0) 870 757 2878	enquiries@berkshire.uk.com
SE Asia	Tel 65 6252 4313	Fax 65 6252 4312	enquiries@berkshire.com.sg
Japan	Tel 81 3 5827 2380	Fax 81 3 5827 2382	master@berkshire.co.jp



#### **Technical Data (In Dry State)**

Attribute		Units	Value	Test Method		
Basis Weight		g/m²	67.6	TAPPI T-410		
Caliper		μm	243	TAPPI T-411		
Fibers	≥100µm	fibers/cm²	200	IEST-RP.CC004.3, Sec 6.1.3 / Sec 6.2.2		
Particles	<u>≥</u> 0.5μm	x10³/cm²	16	IEST-RP.CC004.3, Sec 6.1.3 / Sec 6.2.1		
Sorbency	Capacity	mL/m²	291	IEST-RP.CC004.3, Sec 8.1 modified / Sec 8.2 modified		
	Efficiency	mL/g	4.3			
	Rate	seconds	2			
Non-Volatile Residue	DI Water	g/m²	0.030	IEST-RP.CC004.3, Sec 7.1.2		
	IPA	g/m²	0.0034			
lons	Na+	ppm	49	IEST-RP.CC004.3, Sec 7.2.2		
	K <sup>+</sup>	ppm	4.1			
	Ca <sup>++</sup>	ppm	12			
	Mg <sup>++</sup>	ppm	4.7			
	CI <sup>-</sup>	ppm	30			

#### Notes:

- Technical data represented in this table are typical values at the time of publication. These should not be used as product specifications.
- Due to differences in test methods applied and equipment utilized by different wiper manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions
- Third party testing can be performed upon request

#### **Order Information:**

Product	Number	Size	Shts/pk	Pks/cs	IPA/DI	Saturation	VOC % by Weight	Style
SatPax® 1000	SPX1000.001.12	9x9" (23x23cm)	75	12	70/30	33%	36%	Pouch
SatPax® 1000	SPX1000.002.12	9x9" (23x23cm)	75	12	09/91	45%	6%	Pouch
SatPax® 1000C	SPXC1000.001.12	6x9" (15x23cm)	100	12	70/30	45%	41%	Canister
SatPax® 1000R	SPXC1000.01.12R	6x9" (15x23cm)	100	12	70/30	45%	41%	Refill
SatPax® 1000C	SPXC1000.002.12	6x9" (15x23cm)	100	12	09/91	45%	4.8%	Canister
SatPax® 1000R	SPXC1000.02.12R	6x9" (15x23cm)	100	12	09/91	45%	4.8%	Refill